


LithiX

HERTZ CONTACT INTRAVASCULAR
LITHOTRIPSY CATHETER

LithiX™ Hertz Contact Intravascular Lithotripsy Catheter

TECHNICAL SPECIFICATIONS

LithiX HC-IVL Catheter	
Available Balloon Nominal Diameters	1.5, 2.0, 2.5, 3.0, 3.5 mm
Catheter Working Length	140 cm
Catheter Design	Rapid Exchange
Guiding Catheter Compatibility	6 French (≥ 0.071 " ID or 1.8 mm ID)
Guidewire Compatibility	0.014" or 0.36 mm



Elixir Medical Corporate Headquarters
920 N. McCarthy Blvd, Milpitas, CA 95035 USA
+1.408.636.2000
info@elixirmedical.com

Elixir is a registered trademark and LithiX is a trademark
in the US and internationally for Elixir Medical Corporation.

LithiX Hertz Contact Intravascular Lithotripsy Catheter is CE Mark approved.
©2024 Elixir Medical. All rights reserved. PMN 2050 Rev A

www.ElixirMedical.com

E L I X I R
MEDICAL CORPORATION

LITHiX

HERTZ CONTACT INTRAVASCULAR
LITHOTRIPSY CATHETER

COMPLIANCE CHART

Inflation ATM	Balloon Outer Diameter (mm)				
	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm
5 [NOM]	1.51	2.06	2.57	2.99	3.59
6	1.53	2.09	2.61	3.04	3.65
7	1.55	2.12	2.65	3.09	3.70
8	1.56	2.15	2.69	3.14	3.76
9	1.58	2.18	2.72	3.18	3.80
10	1.59	2.20	2.75	3.21	3.83
11	1.6	2.22	2.77	3.24	3.86
12 [RBP]	1.62	2.24	2.79	3.27	3.88

PRODUCT ORDERING INFORMATION

LithiX Diameter				
1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm
DAR1514	DAR2014	DAR2514	DAR3014	DAR3514

FOR MORE INFORMATION OR TO ORDER THE LITHIX HC-IVL CATHETER:

Contact your local sales representative or email us at Sales@elixirmedical.com

Elixir Medical Corporate Headquarters
920 N. McCarthy Blvd, Milpitas, CA 95035 USA
+1.408.636.2000
info@elixirmedical.com

Elixir is a registered trademark and LithiX is a trademark
in the US and internationally for Elixir Medical Corporation.

LithiX Hertz Contact Intravascular Lithotripsy Catheter is CE Mark approved.
©2024 Elixir Medical. All rights reserved. PMN 2050 Rev A

www.ElixirMedical.com

